


easyPGX[®]

Ready to yo**USE**


diatech
pharmacogenetics



Easy_{PGX}[®] product line - key features

- **Ready to use:** reagents are delivered in 8-well strips preloaded with a complete master mix in a dry, room temperature and stable format
- **Easy to use:** no need for freezing, thawing or pipetting on ice and the few remaining pipetting steps minimize the risk of errors or contaminations
- **High sensitivity:** limit of detection as low as 0.5%
- **Flexible sample requirement:** low DNA input from a variety of sources, including FFPE and plasma
- **Turnaround time:** from tissue to result in less than 3 hours with only 10 minutes of hands-on time
- **Quality assurance:** manufactured under ISO 13485
- **Regulatory:** kits have been designed, developed and validated in accordance with the Directive 98/79/EC on in vitro diagnostic medical devices

From tissue to result in less than 3 hours



System ordering information

Catalog number	Product description	Picture
RT800 CE IVD	Easy _{PGX} qPCR instrument	
RT800-96 CE IVD	Easy _{PGX} qPCR instrument 96	
RT800-SW CE IVD	Easy _{PGX} analysis software	
RT801 CE IVD	Easy _{PGX} dry block	
RT802 CE IVD	Easy _{PGX} centrifuge/vortex 1.5 ml	
RT803 CE IVD	Easy _{PGX} centrifuge/vortex 8-well strips	

Easy_{PGX}[®] ready KRAS
cat. no. RT021 (48 test, CE IVD)

Main features

Detection of the main mutations of exon 2 (codons 12, 13), of exon 3 (codons 59, 61) and of exon 4 (codons 117, 146) of the gene KRAS using 8 oligo mixes. Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene. A specific oligo control mix enables the evaluation of the quality and the quantity of the DNA in each sample.

Controls

The kit is provided with the following controls:

- Positive control sample containing a mixture of synthetic DNA sequences that correspond to each mutation detected by this kit in a background of wild-type genomic DNA.
- Negative control.

Sensitivity

The kit allows the detection of low percentages of mutated allele in presence of high amounts of wild-type genomic DNA by real-time amplification with sequence-specific probes marked with FAM and HEX. **LOD down to 0.5%.**

Starting material

The kit allows the analysis of DNA extracted from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready BRAF
cat. no. RT022 (48 test, CE IVD)

Main features

Detection of the main mutations of codon 600 of the gene BRAF using 4 oligo mixes. Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene. A specific oligo control mix enables the evaluation of the quality and the quantity of the DNA in each sample.

Controls

The kit is provided with the following controls:

- Positive control sample containing a mixture of synthetic DNA sequences that correspond to each mutation detected by this kit in a background of wild-type genomic DNA.
- Negative control.

Sensitivity

The kit allows the detection of low percentages of mutated allele in presence of high amounts of wild-type genomic DNA by real-time amplification with sequence-specific probes marked with FAM and HEX. **LOD down to 0.5%.**

Starting material

The kit allows the analysis of DNA extracted from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.



Easy_{PGX}[®] ready EGFR
cat. no. RT023 (48 test, CE IVD)

Main features

Detection of the main mutations of exons 18, 19, 20, 21 of EGFR gene using 8 oligo mixes. Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene. A specific oligo control mix enables the evaluation of the quality and the quantity of the DNA in each sample.

Controls

The kit is provided with the following controls:

- Positive control sample containing a mixture of synthetic DNA sequences that correspond to each mutation detected by this kit in a background of wild-type genomic DNA.
- Negative control.

Sensitivity

The kit allows the detection of low percentages of mutated allele in presence of high amounts of wild-type genomic DNA by real-time amplification with sequence-specific probes marked with FAM and HEX. **LOD down to 0.5%.**

Starting material

The kit allows the analysis of DNA extracted from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready NRAS
cat. no. RT024 (48 test, CE IVD)

Main features

Detection of the main mutations of exon 2 (codons 12, 13), of exon 3 (codons 59, 61) and of exon 4 (codons 117, 146) of NRAS gene using 8 oligo mixes. Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene. A specific oligo control mix enables the evaluation of the quality and the quantity of the DNA in each sample.

Controls

The kit is provided with the following controls:

- Positive control sample containing a mixture of synthetic DNA sequences that correspond to each mutation detected by this kit in a background of wild-type genomic DNA.
- Negative control.

Sensitivity

The kit allows the detection of low percentages of mutated allele in presence of high amounts of wild-type genomic DNA by real-time amplification with sequence-specific probes marked with FAM and HEX. **LOD down to 0.5%.**

Starting material

The kit allows the analysis of DNA extracted from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.



Easy_{PGX}[®] ready DPYD
cat. no. RT026 (48 test, CE IVD)

Main features

Detection, by allelic discrimination, of the DPYD gene polymorphisms DPYD*2A (IVS14+1G>A, c.1905+1G>A, rs3918290), DPYD*13 (c.1679T>G, rs55886062), DPYD D949V (c.2846A>T, rs67376798) and DPYD IVS10 (c.1129-5923C>G, rs75017182), associated with the toxicity due to the treatment with Fluoropyrimidines, using 4 oligo mixes. Each mix allows the co-amplification of the mutant sequence (FAM) as well as the wild-type sequence (HEX).

Controls

The kit is provided with the following controls:

- DPYD WT positive control: Positive control DNA containing a mixture of synthetic wild-type DNA sequences for the DPYD polymorphisms analyzed.
- DPYD MT positive control: Positive control DNA containing a mixture of synthetic mutant DNA sequences for the DPYD polymorphisms analyzed.
- Negative control.

Starting material

The kit allows the analysis of genomic DNA extracted from whole blood.

Easy_{PGX}[®] ready UGT1A1
cat. no. RT027 (48 test, CE IVD)

Main features

Detection, by allelic discrimination, of the UGT1A1 gene polymorphisms UGT1A1*1 (TA)₆, UGT1A1*28 (TA)₇, UGT1A1*36 (TA)₅ and UGT1A1*37 (TA)₈, associated with the toxicity due to the treatment with Irinotecan, using 1 oligo mix. UGT1A1 mix contains HEX labeled probes for UGT1A1*28 and UGT1A1*37 and FAM labeled probes for UGT1A1*1 and UGT1A1*36.

Controls

The kit is provided with the following controls:

- UGT1A1 WT positive control: Positive control DNA containing synthetic wild-type UGT1A1*1/*1 DNA sequence.
- UGT1A1 MT positive control: Positive control DNA containing a synthetic mutant UGT1A1*28/*28 DNA sequence.
- Negative control.

Starting material

The kit allows the analysis of genomic DNA extracted from whole blood.



Easy^{PGX}[®] ready Thyroid cat. no. RT028 (48 test, CE IVD)

Main features

Detection of the main mutations of exon 2 (codons 12,13), of exon 3 (codons 61) of the genes KRAS, NRAS, HRAS and of the codons 600 and 601 of the gene BRAF using 8 oligo mixes. Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene.

Controls

- The kit is provided with the following controls:
- Positive control sample containing a mixture of synthetic DNA sequences that correspond to each mutation detected by this kit in a background of wild-type genomic DNA.
 - Negative control.

Sensitivity

The kit allows the detection of low percentages of mutated allele in presence of high amounts of wild-type genomic DNA by real-time amplification with sequence-specific probes marked with FAM and HEX. **LOD down to 0.5%.**

Starting material

The kit allows the analysis of DNA extracted from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues, plasma*, and cytological samples.

Helix[®] circulating Nucleic Acid cat. no. H8040 (50 test, CE IVD)

Main features

The kit allows the manual extraction of circulating free DNA (cfDNA) from plasma. The kit Helix[®] Circulating Nucleic Acid, in association with the kit Easy^{PGX}[®] ready EGFR, enables the mutational analysis of EGFR gene in the circulating tumor DNA (liquid biopsy) when the tumor tissue is not evaluable, according to the EMA/129677/2014 recommendations of September 25th 2014. DNA capture by silica membrane and vacuum-based system. The system to concentrate the final eluate up to 3 times is included in the kit.

Starting material

1-5 ml of fresh or frozen plasma.

Turn around time

3 hours

File: Cancers - Tools

Experiment

Name: Administrator

Experiment description:

Add description here:

Platform and application

Select instrument:

RT028 readyPGX ready (CE IVD)

Select kit:

RT028 EasyPGX ready (CE IVD)

Singleplex data and analysis

Select data:

C:\Run\20170808_RT028_ESP1.txt

[OK] [Cancel] [Save]

Analysis of the reaction controls:

Run	Sample	Primer	CA	AM	Result	CA	AM	Result	CA	AM	Result	CA	AM	Result	CA	AM	Result	CA	AM	Result			
1	Positive	FAM	24.0	87.0	OK	28.0	83.0	OK	20.2	22.0	OK	22.8	20.0	OK	21.0	21.0	OK	21.0	21.0	OK	24.0	87.0	OK
2	Negative	HEX	24.0	87.0	OK	28.0	83.0	OK	20.2	22.0	OK	22.8	20.0	OK	21.0	21.0	OK	21.0	21.0	OK	24.0	87.0	OK

Analysis of the sample control mix and mutation assay:

Run	Sample	Primer	CA	AM	Result	CA	AM	Result	CA	AM	Result	CA	AM	Result	CA	AM	Result	CA	AM	Result			
1	GGFR_FFPE_sample 1	FAM	4	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK
2	GGFR_FFPE_sample 2	HEX	22	220	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK
3	GGFR_FFPE_sample 3	FAM	4	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK
4	GGFR_FFPE_sample 4	HEX	22	220	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK
5	GGFR_FFPE_sample 5	FAM	4	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK
6	GGFR_FFPE_sample 6	HEX	22	220	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK
7	GGFR_FFPE_sample 7	FAM	4	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK
8	GGFR_FFPE_sample 8	HEX	22	220	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK
9	GGFR_FFPE_sample 9	FAM	4	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK	11	WT	OK
10	GGFR_FFPE_sample 10	HEX	22	220	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK	12	86	OK

Warnings

Code	Description
001	Invalid error in the set up of the reaction run. It is not possible to analyze the samples (see Troubleshooting).
002	Invalid concentration if it is not possible to analyze the samples (see Troubleshooting).
003	Unoptimal amount of starting DNA on PCR inhibition (see Troubleshooting).
004	Excess of DNA. Sample must be diluted with water so that Ca be in the range indicated (see Troubleshooting).
005	Not sufficient template / PCR inhibition / inhibitor during sample's denaturation (see Troubleshooting).
006	Variable amount of DNA. Proceed with the analysis of the sample.

easyPGX

* Please note that extraction from plasma is sold separately (cat.n. H8040)

For information please contact:

diatech pharmacogenetics

Diatech Pharmacogenetics srl

Via Ignazio Silone 1b - 60035 Jesi (An) Italy

Phone +39 0731 213 243

export@diatechpgx.com

www.diatechpharmacogenetics.com